

Full-element Analysis

Rapid Accurate Non-destructive


Skyray
Skyray Instrument



Cement




Steel and other
nonferrous
metals



RoHS



Aggregate



EDX 3600B

X-ray Fluorescence
Spectrometer

Skyray Instrument

❖ **USA Office:**
Add: 639 Granite Street, suite 222 Braintree, MA 02184
Tel: 617-202-3879
Fax: 617-202-3878
Website: www.skyrayinstrument.com
E-mail: sales@skyray-instrument.com

❖ **Kunshan Office:**
Add: Skyray Building, Tsinghua Science Park, 1666,
South Weicheng Rd., Kunshan, Jiangsu Province
Tel: 0512-57017006 57017018 57017886 57017880
Fax: 0512-57017357
Website: www.gold-tester.com
E-Mail: sales@gold-tester.com

❖ **Shenzhen Office:**
Add: Floor 3(rd), Block D, Keepon Industry Zone, Chongqing
Rd., Fuyong Town, Bao'an District, Shenzhen, China
Tel: 0755-81459355 81459357 81459359
Fax: 0755-81459301
Hotline: 800-9993-800


X-ray Fluorescence Spectrometer

Skyray Skyray Instrument

Full-element Analysis Expert

EDX 3600B

Engaged in cement, steel, nonferrous metals, aggregate and RoHS detection



Rapid
Accurate
Non-destructive

The most effective world-leading X-ray generator and high-voltage unit make the instrument produce results compliant with State standards.

The instrument is specially designed for metallurgy and building material industry.

Precision Instruments
Made in Skyray

XRF technology is used in EDX3600B X-ray Fluorescence spectrometer for rapid and accurate elemental analysis. The technology features low-energy X-ray which generates good excitation results of light elements such as Si, S, Na and Mg. Moreover, with short test time, the test efficiency has been significantly improved. EDX3600B has good energy linearity, energy resolution, spectrum property and high peak-background ratio by using UHRD detector. It is consistent due to the automatic spectrum stabilizing device. The original spectrum can be easily decoupled by applying UHRD technology, which improves the measured analytical precision of light elements of Si, S, Al, etc. Owing to multi-parameter linear regression method, the absorption and repelling effects between elements can be highly reduced.



Performance traits

- It performs professional full-element analysis on cement, steel, aggregate, plating thickness detection and RoHS.
- In-built SNE improves the signal processing ability up to 25 times.
- The collimators and filters can be switched automatically for different samples.
- Electro-cooling UHRD detector instead of liquid nitrogen cooling detector
- Intelligent full-element analysis software matches with the hardware well.

Product Parameters


Name: Skyray X-ray Fluorescence Spectrometer
Model: EDX3600B
Input Voltage: AC 110V/220V
Consumed Power: 200W
Ambient Temperature: 15°C-30°C
Ambient Humidity: 35%-70%
Area of Sample Chamber: 320mm×180mm
Size: 650mm×608mm×466mm
Weight: 75kg

Unique configurations

- Signal-to-Noise Enhancer (SNE)
- Light path enhancement system
- Electro-cooling UHRD detector
- In-built high resolution camera
- Automatic collimator and filter switch
- Precise mobile platform
- Enhanced metal sensitivity analyzer

Technical specifications

Range of measurable elements: Na to U
Range of element content: 1ppm-99.99%
Ability of simultaneous analysis: 24 elements
Functions: full-element analysis of cement, steel and aggregate
Plating thickness: more than 11 layers, up to 0.005um each layer
Analysis accuracy: 0.05%
Forms of samples: powder, solid and liquid
Testing time: 60-200s
Energy resolution: (150±5) eV
Tube voltage: 5kV-50kV
Tube current: 50uA-1000uA



EDX 3600B

EDX 3600B



Vacuum system- a technological advance

Compared with tradition instruments, this one is more advantageous. It shields against the atmospheric influence and greatly enlarges the measurement scope. By a touch of the finger, the instrument opens and shuts, rises and falls.

The inner structure- a perfect combination of strength and beauty.

The whole steel structure is a guarantee of the strength reliability, while the external appearance owns its smoothness and elegance to plastic material.

Automatic collimator and filter switching- a tradition of intelligence system

This exempts you from the complicated manual operation and its technological charm blooms in the humanity consideration.

Triple safety protection mode

Independent matrix effect correction models

Multi-variable non-linear regression procedure

Arbitrary optional analysis and identification models

X-ray Fluorescence Spectrometer
EDX 3600B

Furthermore, EDX3600B is used in metallurgical and mineral industries. Our customers are:

NO.	Name of the organization	Application scope
1	National Museum of China	To test ceramics, ancient pottery, bronze ware or help identify the age and validate their authenticity.
2	Shanxi Heavy Vehicle Group Co., Ltd.	To test iron and steel, screws, alloys and so on. To judge whether they are up to the standard.
3	Zhejiang Zhengtai Electrical Appliance Co., Ltd.	To test alloys, steel, stainless steel, plating things and other ironware.
4	Henan Baige Group Co., Ltd.	To test the components of corundum finished products and control the quality of the production process.
5	Jiangmen Fuyi Magnetic Material Co., Ltd.	To test the mixture ratio of ferrite materials and magnetic materials.
6	Jiangsu Lianyungang Geography Engineering Exploration Institute	To test Niobium and Tantalum ore. or used for mineral inspection and mineral collection.
7	Shandong University of Science and Technology.	Used for teaching and scientific study, for example, minerals and alloys test.

The difference between the application of X-ray analytical instrument EDX3600B and traditional chemical method in iron and steel industry:

The iron and steel industry is a manufacturing industry running on large scale and ceaselessness. Raw materials, accessories and middle-products are required to be quickly tested before the results is used for guiding the production process. For such a characterized industry, rapidness is a must in getting the test result. So the analysis of this industry has two distinct traits: samples are various and the tests are in high speed.

Item	Traditional chemical method	EDX3600 family
Speed	Slow, 10 ~ 30 min.	Quick, 2 ~ 5min.
Analytical result	Greatly influenced by operator, poor repeatability	No operator influence, high precision and repeatability
Labor intensity	High	Low, most automatically done
Elements analyzed	Only one element at a time	Tens of elements at a time
Chemical relativity	Analysis method varies with different sample chemical property	Pure physical method, no relativity with sample chemical property
Cost	High	Low
Skill requirement	High, long-time training	Low, simple training

Applications

EDX3600B X-ray spectrometer has been widely used in the iron and steel industry and warmly welcomed by our customers.


Some of our customers are:

NO.	Name of the organization	Application
1	Shanxin Hanzhong Steel Co., Ltd.	To test blast furnace slag, pig iron, sinter, pellet, dolomite, bentonite, etc.
2	Xiaojiao Steel Rolling Plant, Fujian Sanming Steel Group Co., Ltd.	To test pig iron, sinter, limestone, blast furnace slag, etc.
3	Daguanshan Mine, Taiyuan Iron & Steel (Group) Co., Ltd.	To test quartz, dolomite, ordinary silicon and special silicon, etc.
4	Tangshan Baoye Group Iron and Steel Co., Ltd.	To test pig iron, blast furnace slag, sinter, iron concentrate fines, fused Magnesia, etc.

EDX 3600B


Examples

Cement industry



Element	content	intensity
Na	0.321	01.610
Mg	2.954	99.550
Al	4.042	182.950
Si	14.939	750.010
K	0.593	326.660
Ca	41.094	8114.300
Fe	0.213	62.100

Working Principle



Characteristic X-radiation of element

Each element will emit X-ray at its own energy level when excited. This X-ray is characteristic and called X-ray fluorescence. It is the foundation of analysis.

Scattering

It is the background of spectrum.

Photoelectron

The photoelectron is the foundation of detector. In the sample, the X-ray intensity of every element is expressed as I₁, I₂, I₃, I₄, I₅ ... respectively. The element content C is the function of X-ray fluorescence intensity I, expressed as following:


$$C = f(I_1, I_2, I_3, I_4, I_5, \dots)$$

This equation is too complicated and can be simplified as:

$$C = K_1 I_1 + K_2 I_2 + K_3 I_3 + K_4 I_4 + K_5 I_5 \dots$$

Where

C is the element content in the sample; I₁, I₂, I₃, I₄, I₅ ... are X-ray intensity of element respectively; K₁, K₂, K₃, K₄, K₅ ... are coefficients which can be determined by measuring known standard sample to calibrate.



Analytical data of light elements taken by EDX3600B

Times	Elements Results	Ca content	Na content	Mg content	Al content	Si content	K content	Ti content	Fe content
Test1		44.410	0.458	3.563	4.582	10.012	0.344	0.113	3.388
Test2		44.400	0.448	3.460	4.442	9.801	0.352	0.113	3.385
Test3		44.448	0.427	3.480	4.493	9.954	0.347	0.112	3.368
Test4		44.461	0.430	3.652	4.397	9.804	0.349	0.111	3.389
Test5		44.423	0.448	3.608	4.491	9.879	0.342	0.113	3.388
Test6		44.475	0.450	3.534	4.519	9.917	0.347	0.112	3.375
Test7		44.446	0.444	3.782	4.663	9.934	0.351	0.111	3.389
Test8		44.347	0.446	3.729	4.442	10.110	0.350	0.109	3.365
Test9		44.355	0.450	3.560	4.450	10.090	0.352	0.111	3.354
Test10		44.444	0.452	3.581	4.359	10.076	0.350	0.113	3.374
Test11		44.444	0.464	3.554	4.548	9.688	0.351	0.113	3.372
Test12		44.429	0.448	3.655	4.499	9.789	0.352	0.113	3.350
Test13		44.430	0.453	3.717	4.521	10.093	0.350	0.113	3.359
Test14		44.490	0.451	3.779	4.466	10.206	0.348	0.111	3.362
Test15		44.462	0.446	3.679	4.483	10.212	0.348	0.115	3.362
Test16		44.494	0.445	3.445	4.426	10.044	0.355	0.109	3.350
Test17		44.440	0.435	3.720	4.418	10.053	0.348	0.111	3.355
Test18		44.418	0.438	3.640	4.543	10.067	0.354	0.108	3.369
Test19		44.428	0.427	3.705	4.359	10.310	0.350	0.112	3.371
Test20		44.551	0.466	3.646	4.659	10.391	0.350	0.110	3.368
Average value of content		44.44	0.45	3.62	4.49	10.02	0.35	0.11	3.37
Standard deviation of the measurement S _n		0.04	0.01	0.10	0.08	0.18	0.00	0.00	0.01
3s value		0.13	0.03	0.29	0.25	0.53	0.01	0.01	0.04
Relative standard deviation RSD (%)		0.10%	2.34%	2.71%	1.82%	1.76%	0.86%	1.51%	0.38%

$$S_{(n)} = \sqrt{\frac{\sum_{i=1}^n (N_i - \bar{N})^2}{n-1}}$$

$$\bar{N} = \frac{\sum_{i=1}^n N_i}{n}$$

S_n --- Standard deviation after n times of tests
N --- times of test

\bar{N} --- Average value after n times of tests
RSD can be calculated by using the equation below

$$RSD = \frac{S_{(n)}}{\bar{N}} \times 100\%$$

Sample name: Cement Xs2 #

EDX3600B


Iron and steel industry

Pig iron

Pig iron is the most important middle product in the production process of iron and steel industry. Its quality is directly concerned with the quality of the final product. It determines the usage amount of energy and accessories before it is processed to steel billet. Si, Mn, P and S in the pig iron can also be detected by X-ray fluorescence spectrometer. The test is highly precise.


The spectrum and working curve of real iron test taken by EDX3600B are listed below.

The pig iron spectrum taken by EDX3600B



The spectrum figure shows that the spectra line of Si, Mn, P and S are very clear. The sample element measurement can be obtained easily from the spectrum. Moreover, the spectra line of the accessory elements such as Al and Cr are also very clear. They can also be detected at the same time.

Working curve and relative data of Pig iron-Si




Spectrum name	Content	Intensity
LQ-6	1.16	172.06
13662Q	0.39	89.08
11049Q	2.37	348.8
11051Q	2.45	372.8
11048Q	2.19	337.96
11049Q	2.52	391.27
4265-1	0.25	66.77
4265-2	0.28	73.63
14297	0.46	116.53
14299-1	0.47	82.44
14300-1	0.55	104.26
14300-2	0.69	122.17
14302-2	0.36	76.36
14323-1	0.37	94.36
14361-1	0.33	73.13

The test accuracy of the actual pig iron test

Elements	Si	Mn	P	S
Range of content	0.1-1.0	0.1-1.0	0.05-0.1	0.010-0.10
(SD) Test accuracy	0.02	0.015	0.004	0.002

Sintering ore




During the test of sintering ore, EDX3600B mainly detects the elements like Fe, Ca and Si. The test accuracy is very reliable. It can also detect the accessory elements like Mg, P, S, Cr, Mn and Ti.


The spectrum and working curve of real sintering ore taken by EDX3600B are listed below.

Sintering ore spectrum taken by EDX3600B

A. spectrum of high-alkali sintering ore




B. Spectrum of low-alkali sintering ore



The spectrum figure shows that the spectra line of Fe, Ca and Si are very clear. The sample element measurement can be obtained easily from the spectrum. Moreover, the spectra line of the accessory elements such as Al, K, S, Fe and Mg are also very clear. They can also be detected at the same time.

The working curve of Ca in sintering ore and other relative data



Spectrum name	Content	Intensity
YSBC15704	15.79	3838.38
CS-D8	2.09	647.39
LS-24	13.13	3344.89
LS-22	12.48	3186.29
LSD-10	2.28	661.17
LSD-11	2.46	712.63
S4-23	12.77	3264.87
LS-25	12.99	3283.42
CSD-9	2.16	648.77
LS-23	12.69	3216.76
LS-21	12.84	3296.33
GBW07219A	11.03	2817.66
S4-45	13.38	3366.53

The test accuracy of the actual test of sintering ore

Elements	Ca	Si	Fe
Range of content	5.00-20.00	5.01-15.00	≥50.00
(SD) Test accuracy	≤0.15	≤0.15	≤0.2



Full-element Analysis Expert

X-ray Fluorescence Spectrometer

Skyray Skyray Instrument


EDX3600B

Ferrosilicon


Ferrosilicon is usually used to test the content of Si, Fe, and Ca. On our instrument, it can test the elements like Al, Ca, Mn, Cr, and Ti with great precision.

Below are the spectrum and working curves in a test of ferrosilicon by EDX3600B

The spectrum of ferrosilicon taken by EDX3600B



The working curve and relative data of Mn in ferrosilicon



Spectrum name	Content	Intensity
BH0225-1	0.411	150.21
BH0225-2	0.308	121.09
BH0301-4	0.17	102.89
YSBC28604	0.131	96.09
YSBC28607	0.13	96.51
第2001号	0.26	110.34
YSBC28608	0.13	94.36
第2002号	0.26	113.24
BH1913-1	0.128	93.97
BH0301-6	0.14	93.23
GBW01422A	0.172	97.64
YSB14602	0.048	72.52

The accuracy of the actual test of ferrosilicon

Elements	Si	Ca	Fe	Mn
Range of content	≥ 50.00	0.10-2.00	20.00 - 50.00	≤ 0.5
(SD) Test accuracy	≤ 0.2	≤ 0.015	≤ 0.10	0.01

PERIODIC TABLE OF ELEMENTS (Characteristic X-Ray energy table)

Skyray Skyray Instrument

Hotline: 800-9993-800

Address: Skyray Building, Tsinghua Science Park, 1666, South Weicheng Rd., Kunshan, Jiangsu Province
Tel: 0512-57017006 57017018 57017886 57017880

IA	IIA	IIIB	IVB	VB	VIB	VII B	VIII	IB	II B	0							
1 H 1.008										2 He 4.008							
3 Li 6.94 0.052	4 Be 9.012 0.110																
11 Na 22.99 1.041	12 Mg 24.31 1.254																
19 K 39.1 3.312 3.589	20 Ca 40.08 3.690 0.341	21 Sc 44.96 4.088 0.395	22 Ti 47.90 4.508 0.452	23 V 50.94 4.949 0.51	24 Cr 51.99 5.411 0.571	25 Mn 54.94 5.895 0.636	26 Fe 55.84 6.400 0.704	27 Co 58.93 6.925 0.775	28 Ni 58.7 7.472 0.849	29 Cu 63.54 8.041 0.928	30 Zn 65.38 8.631 1.009	31 Ga 69.72 9.243 1.096	32 Ge 72.5 9.876 1.186	33 As 74.92 10.66 1.261	34 Se 78.96 11.223 1.375	35 Br 79.90 11.91 1.48	36 Kr 83.8 12.63 1.587
37 Rb 85.47 13.38 14.97 1.694 1.752	38 Sr 87.82 14.14 15.85 1.806 1.872	39 Y 88.91 14.93 16.75 1.922 1.996	40 Zr 91.22 15.75 17.69 2.042 2.124	41 Nb 92.91 16.58 18.65 2.166 2.257	42 Mo 95.94 17.44 19.63 2.293 2.395	43 Tc #(99) 18.33 20.65 2.424 2.538	44 Ru 101.0 19.24 21.69 2.558 2.683	45 Rh 102.9 20.17 22.76 2.696 2.834	46 Pd 106.4 21.12 23.86 2.838 2.990	47 Ag 107.9 22.10 24.99 2.984 3.151	48 Cd 112.4 23.11 26.14 3.133 3.316	49 In 114.8 24.14 27.38 3.287 3.487	50 Sn 118.6 25.19 28.60 3.444 3.662	51 Sb 121.7 26.27 29.85 3.605 3.843	52 Te 127.6 27.8 31.75 3.789 4.028	53 I 126.9 28.51 32.44 3.937 4.220	54 Xe 131.3 29.67 33.78 4.111 4.422
55 Cs 132.9 30.85 35.15 4.286 4.620 5.280 3.794	56 Ba 137.3 32.07 36.55 4.467 4.828 5.531 3.953	Ln	72 Hf 178.4 55.38 63.56 7.898	73 Ta 180.9 57.11 65.56 8.145	74 W 183.8 58.86 67.59 8.396	75 Re 186.2 60.66 69.66 8.651	76 Os 190.2 62.48 71.78 8.910	77 Ir 192.2 64.35 73.93 9.173	78 Pt 195.0 66.25 76.13 9.441	79 Au 197.0 68.19 78.37 9.711	80 Hg 200.5 70.16 80.66 9.987	81 Tl 204.3 72.18 82.99 10.27	82 Pb 207.2 74.23 85.36 10.55	83 Bi 208.9 76.32 87.77 10.84	84 Po #(209) 78.46 90.24 11.13	85 At #(210) 80.64 92.75 11.42	86 Rn #(222) 82.86 95.32 11.72
87 Fr #(223) 82.12 97.93 12.03 14.77	88 Ra 226.0 87.44 100.6 12.34 17.8 10.60	An	72 Hf 178.4 55.38 63.56 7.898	73 Ta 180.9 57.11 65.56 8.145	74 W 183.8 58.86 67.59 8.396	75 Re 186.2 60.66 69.66 8.651	76 Os 190.2 62.48 71.78 8.910	77 Ir 192.2 64.35 73.93 9.173	78 Pt 195.0 66.25 76.13 9.441	79 Au 197.0 68.19 78.37 9.711	80 Hg 200.5 70.16 80.66 9.987	81 Tl 204.3 72.18 82.99 10.27	82 Pb 207.2 74.23 85.36 10.55	83 Bi 208.9 76.32 87.77 10.84	84 Po #(209) 78.46 90.24 11.13	85 At #(210) 80.64 92.75 11.42	86 Rn #(222) 82.86 95.32 11.72

■ Alkali Metals ■ Non-Metal ■ Transitional element
■ Halogen ■ Lanthanoids ■ Actinides
■ Alkaline earth ■ Rare gases ■ Main group metal

Mark 1: #Radioactive Elements *Man Made Elements
2: All the Numbers are ordered one by one in this way, Atomic Number
Element Symbol, Atomic Weight, K α , K β , L α , L β , L γ , L ϵ

Ln	57 La	58 Ce	59 Pr	60 Nd	61 Pm	62 Sm	63 Eu	64 Gd	65 Tb	66 Dy	67 Ho	68 Er	69 Tm	70 Yb	71 Lu
	138.9	140.1	140.9	144.2	#(147)	150.4	152.0	157.2	158.9	162.5	164.9	167.2	168.9	173.0	175.0
	33.30	34.57	35.86	37.19	38.54	39.91	41.32	42.76	44.23	45.73	47.26	48.82	50.41	52.04	53.59
	37.99	39.45	40.95	42.48	44.05	45.65	47.28	48.95	50.65	52.38	54.16	55.96	57.81	59.69	61.61
	4.651	4.840	5.034	5.230	5.431	5.636	5.846	6.059	6.275	6.495	6.720	6.948	7.181	7.414	7.654
	5.043	5.262	5.489	5.722	5.956	6.206	6.456	6.714	6.979	7.249	7.528	7.810	8.103	8.401	8.708
	5.789	6.052	6.322	6.602	6.891	7.180	7.478	7.778	8.104	8.418	8.748	9.089	9.424	9.779	10.14
	4.124	4.287	4.452	4.632	4.816	4.994	5.176	5.361	5.546	5.742	5.942	6.152	6.341	6.544	6.752
An	89 Ac	90 Th	91 Pa	92 U	93 Np	94 Pu	95 Am	96 Cm	97 Bk	98 Cf	99 Es	100 Fm			
	#(227)	#(232)	231.0	238.0	(237)	(244)	(243)	(247)	(247)	(251)	(252)	(257)			
	89.79	92.19	94.64	97.14	99.69	102.3	104.9	107.7	110.5	113.3	116.2	119.2			
	103.3	106.1	108.9	111.8	114.7	117.7	120.8	123.9	127.1	130.4	133.7	137.2			
	12.65	12.97	13.29	13.61	13.95	14.28	14.62	14.96	15.31	15.66	16.02	16.38			
	15.71	16.2	16.7	17.22	17.74	18.28	18.83	19.39	19.97	20.56	21.17	21.79			
	18.41	18.98	19.55	20.16	20.77	21.40	22.04	22.69	23.37	24.06	24.76	25.47			
		11.12	11.36	11.62	11.89	12.12	12.38								